


NGS QC

An introduction




“I think you’ll find that mine is bigger...”

Outline



The Big Picture



Common Sequence Formats

FASTA

- Simple
- Nucleotide or amino acid strings
- No quality info
- Compressible (.gz)

FASTQ

- Mildly complex
- Nucleotide strings (not AA)
- Quality information included
- Compressible (.gz)

FAST5

- Complex (HDF5)
- Nanopore Data
- Nucleotide strings (not AA)
- Raw *squiggles*
- Natively compressed

Common Sequence Formats

FASTA

- Simple
- Nucleotide or amino acid strings
- No quality info
- Compressible (.gz)

FASTQ

- Mildly complex
- Nucleotide strings (not AA)
- Quality information included
- Compressible (.gz)

FAST5

- Complex (HDF5)
- Nanopore Data
- Nucleotide strings (not AA)
- Raw *squiggles*
- Natively compressed

The FASTA format

>sequence 1

CATCGATCGCATGCTACTGACTG
CATGCTCGCGCCCCCCCCCGATG

>sequence 2

ACTGACTCGCGCGCGCGGGGG
GAGCTGATGTG

>sequence 3

CATCGATCGCATGCTACTGACTG
CATGCTCGCGCCCCCCCCCGATG
ACTGACTCGCGCGCGCGGGGG
GAGCTGATGTG

The FASTA format

>sequence 1

```
CATCGATCGCATGCTACTGACTG  
CATGCTCGCGCCCCCCCCCGATG
```

“interleaved”

>sequence 2

```
ACTGACTCGCGCGCGCGGGGG  
GAGCTGATGTG
```

>sequence 3

```
CATCGATCGCATGCTACTGACTG  
CATGCTCGCGCCCCCCCCCGATG  
ACTGACTCGCGCGCGCGGGGG  
GAGCTGATGTG
```

The FASTA format

>sequence 1

CATCGATCGCATGCTACTGACTG

CATGCTCGCGCCCCCCCCCGATG

>sequence 2

ACTGACTCGCGCGCGCGGGGGG

GAGCTGATGTG

>sequence 3

CATCGATCGCATGCTACTGACTG

CATGCTCGCGCCCCCCCCCGATG

ACTGACTCGCGCGCGCGGGGGG

GAGCTGATGTG

The FASTA format

“non-interleaved”

>sequence 1

CATCGATCGCATGCTACTGACTGCATGCTCGCGCCCCCCCCCGATG.....

>sequence 2

ACTGACTCGCGCGCGCGGGGGGAGCTGATGTG

>sequence 3

CATCGATCGCATGCTACTGACTGCATGCTCGCGCCCCCCCCCGATGAC...

The FASTQ format

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

The FASTQ format

Sequence ID

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

The FASTQ format

Sequence

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

The FASTQ format

+ description (or empty)

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

The FASTQ format

+ description (or empty)

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

The FASTQ format

Quality score of each base

```
@ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
GTAGAACTGGTACGGACAAGGGGAATCTGACTGTAG  
+ILLUMINA-F6C19_0048_FC:5:1:12440:1460#0/1  
hhhhhhhhhhhhghhhhhhhhehhhedhhhhfhhhhhh
```

Illumina Sequence ID Lines: A Decoder

@M01137:30:000000000-AA299:1:1101:10929:1966

M01137	the unique instrument name
30	the run id
000000000-AA29	the flowcell id
1	flowcell lane
1101	tile number within the flowcell lane
10929	'x'-coordinate of the cluster within the tile
1966	'y'-coordinate of the cluster within the tile
1 or 2 (not shown, optional)	the member of a pair, 1 or 2 (paired-end or mate-pair reads only)
ATCACG (not shown, optional)	index sequence

Quality Scores

- Phred Score
- $Q = -10 \cdot \log_{10} P$ P = probability the base call is incorrect
- ASCII (character) - 33

Phred Quality Score	Probability of incorrect base call	Base call accuracy
0	1	0 %
10	1 in 10	90 %
20	1 in 100	99 %
30	1 in 1000	99.9 %
40	1 in 10000	99.99 %
50	1 in 100000	99.999 %
93	1 in 20000000000	99.9999995 %

Why QC NGS Data?

OPEN ACCESS Freely available online

PLOS ONE

An Extensive Evaluation of Read Trimming Effects on Illumina NGS Data Analysis

Cristian Del Fabbro¹, Simone Scalabrin², Michele Morgante¹, Federico M. Giorgi^{1,3*}

“Trimming is shown to increase the **quality** and **reliability** of the analysis, with concurrent gains in terms of **execution time** and **computational resources** needed”

Types of Trimming

Quality

- remove low quality bases and reads
 - **Q20** (1% error) and **Q30** (0.1% error) are standard
- Remove too short reads
- Too many 'N' (uncalled bases)

Complexity

- simple repeats (e.g. TGTGTGTGTGTG)
- Homopolymers (e.g., AAAAAAAAAA)


Contamination

- Sequencing adapters!!!!!!
- lab contamination (human, bacteria)
- Environmental contamination


Low Quality Sequences Before Trimming (Puma 454 sequences)





Same Sequences After Trimming (Puma 454 sequences)




Types of Trimming




Types of Trimming




Types of Trimming




Sliding Window Trimming




Types of Trimming




Sliding Window Trimming




Types of Trimming




Sliding Window Trimming




Types of Trimming



Sliding Window Trimming



Adapter Contamination



Adapter Contamination

❌ Overrepresented sequences

Sequence	Count	Percentage	Possible Source
GACTAAGCAGTGGTATCAACGCAGAGTACATGGGGACACTTGTTTCTGAC	19391	5.415739186535921	No Hit
GACTAAGCAGTGGTATCAACGCAGAGTACATGGGGACACTTGCTTCTGAC	11325	3.162974900083508	No Hit
GACTAAGCAGTGGTATCAACGCAGAGTACATGGGACACTTGTTTCTGACA	9229	2.5775801636088915	No Hit
GACTAAGCAGTGGTATCAACGCAGAGTACATGGGGACACTTGTTTCTGAC	6443	1.7881850822373475	No Hit

Download [Graphics](#)

gn||uv|NGB00593.1:1-30 Evrogen Mint PlugOligo-1 adapter

Sequence ID: Length: 30 Number of Matches: 1

Range 1: 1 to 25 [Graphics](#)

▼ Next Match ▲ Previous Match

Score	Expect	Identities	Gaps	Strand
50.6 bits(25)	0.001	25/25(100%)	0/25(0%)	Plus/Plus

Query	5	AAGCAGTGGTATCAACGCAGAGTAC	29
Sbjct	1	AAGCAGTGGTATCAACGCAGAGTAC	25

Error Correction (Illumina data)

GAGE: A critical evaluation of genome assemblies and assembly algorithms

Steven L. Salzberg,^{1,7} Adam M. Phillippy,² Aleksey Zimin,³ Daniela Puiu,¹ Tanja Magoc,¹ Sergey Koren,^{2,4} Todd J. Treangen,¹ Michael C. Schatz,⁵ Arthur L. Delcher,⁶ Michael Roberts,³ Guillaume Marçais,³ Mihai Pop,⁴ and James A. Yorke³

“For all four genomes and for all eight assemblers used in GAGE, the best assemblies were created from reads that had been processed through extensive error correction routines”

Illumina Sequencing Errors: ~0.1 - 1%, Substitution errors

Error Correction



High quality base, but an Error

Error Correction: *K*-mer Counting

$k = 4$ AGCTGTGG

Error Correction: *K*-mer Counting

$k = 4$

AGCTGTGG



AGCT

Error Correction: *K*-mer Counting

$k = 4$

AGCTGTGG



AGCT

GCTG

Error Correction: *K*-mer Counting

$k = 4$

AGCTGTGG



AGCT

GCTG

CTGT

Error Correction: *K*-mer Counting

$k = 4$

AGCTGTGG



AGCT

GCTG

CTGT

TGTG

Error Correction: *K*-mer Counting

$k = 4$

AGCTGTGG



AGCT

GCTG

CTGT

TGTG

GTGG

Error Correction: *K*-mer Counting

$k = 6$ AGCTGTGG

Error Correction: *K*-mer Counting

k = 6

AGCTGTGG



AGCTGT

Error Correction: *K*-mer Counting

$k = 6$

AGCTGTGG



AGCTGT
GCTGTG

Error Correction: *K*-mer Counting

$k = 6$

AGCTGTGG




AGCTGT

GCTGTG

CTGTGG

Error Correction: *K*-mer Counting


- Expected Distribution of *k*-mer frequency



DSK; Rizk et al. 2013

Error Correction: *K*-mer Counting

- Expected Distribution of *k*-mer frequency




DSK; Rizk et al. 2013

Error Correction: *K*-mer Counting

- Expected Distribution of *k*-mer frequency

Corrected




DSK; Rizk et al. 2013

Error Correction: *K*-mer Counting

- Expected Distribution of *k*-mer frequency

Estimate
genome
size




DSK; Rizk et al. 2013

Error Correction: *K*-mer Counting

- Expected Distribution of *k*-mer frequency

Estimate
genome
size




DSK; Rizk et al. 2013

$G = C / P$
 G = genome size
 C = total count of
true *k*-mers
 P = peak coverage

Error Correction: *K*-mer Counting


- Expected Distribution of *k*-mer frequency

Estimate
repetitive
content



DSK; Rizk et al. 2013

K-mer Profiling Example: *Raillietiella orientalis*



Recap: NGS QC

Remove low quality bases and reads

Identify and remove adapter contamination

Optional: Correct substitution sequencing errors

Optional: De-duplication

To Your Terminals!

